



**KENT
SURGICAL**

Baskets, Stringers & Accessories

Version : 1




Scan this QR code to explore our website and learn more about our products and services.
We invite you to discover how we can meet your needs and exceed your expectations.

Kent Surgical Sterilization Baskets: Secure and Efficient Instrument Handling

Kent Surgical, alongside its extensive range of surgical instruments, offers a comprehensive selection of sterilization baskets designed to safely hold instruments during cleaning, sterilization, and transport. These baskets are manufactured from high-quality AISI 304 stainless steel and boast a superior electropolished burr-free finish, ensuring optimal hygiene and protection for your instruments.


Key Features:

- **Durable and Long-lasting:** Constructed from robust AISI 304 stainless steel for enduring use.
- **Superior Finish:** Electropolished finish eliminates burrs, minimizing instrument damage and facilitating cleaning.
- **Compatible with Various Sterilization Methods:** Suitable for both steam and gas plasma sterilization.



Perforated Sheet Flat Base (FB)

FB Baskets are ideal for instrument sterilization with side handles.
Mesh Size Approx. 4 mm x 4 mm (bottom) / 2 mm x 2 mm (side)



FEET FOR FB TRAY

	Size in mm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
A	480 x 240 x 70	18.90 x 9.45 x 2.76	FB482407
	480 x 250 x 50	18.90 x 9.84 x 1.97	FB482405
B	400 x 240 x 70	15.75 x 9.45 x 2.76	FB402407
	400 x 240 x 50	15.75 x 9.45 x 1.97	FB402405
C	240 x 240 x 70	9.45 x 9.45 x 2.76	FB242407
	240 x 240 x 50	9.45 x 9.45 x 1.97	FB242405
	Feet for FB Tray	Adds an additional 13mm in total height	FB Feet

Multipurpose Wire Mesh (WM)

WM Baskets are ideal for cleaning and sterilization with side handles.
Mesh size Approx. 4 mm x 4 mm



A



B



C





	Size in mm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
A	480 x 250 x 70	18.90 x 9.84 x 2.76	WM482507
	480 x 250 x 50	18.90 x 9.84 x 1.97	WM482505
B	405 x 250 x 50	15.94 x 9.84 x 1.97	WM402505
C	240 x 250 x 70	9.45 x 9.84 x 2.76	WM242507
	240 x 250 x 50	9.45 x 9.84 x 1.97	WM242505

Micro Mesh Baskets (MM)

MM Baskets are ideal for small instruments they come with Lid for added protection.

Mesh Size Approx. 1.6 mm x 1.6 mm



	Size in mm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
A	90 x 90 x 50	3.54 x 3.54 x 1.97	MM090905
B	160 x 40 x 20	6.30 x 1.57 x 0.79	MM160402
C	225 x 110 x 50	8.86 x 4.33 x 1.97	MM221105

Ultra Mesh Basket (UM)

UM Basket are ideal for very fine instruments they come with Lid for added Protection.


Mesh Size Approx. 0.5 mm x 0.5 mm



	Size in mm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
A	40 x 40 x 20	1.57 x 1.57 x 0.79	UM040402
B	80 x 40 x 20	3.15 x 1.57 x 0.79	UM080402
C	80 x 80 x 40	3.15 x 3.15 x 1.57	UM080804

Endoscopy Basket (EB)

Endoscopy Baskets are specifically designed for rigid endoscopes featuring a fully removable lid, silicon inserts, and a locking mechanism for secure storage and transport.



	Size in mm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
A	180 x 80 x 52	7.09 x 3.15 x 2.05	EB180805
B	290 x 80 x 52	11.42 x 3.15 x 2.05	EB290805
C	460 x 80 x 52	18.11 x 3.15 x 2.05	EB460805
D	640 x 150 x 77 (With Micro-Mesh Basket)	25.20 x 5.91 x 3.03	EB641508

Stringers

Kent Surgical Stringers are made from high-grade stainless steel and offers the perfect solution for keeping your instruments organized and protected. Three Models to choose from PL, PS, and US.



Pull-to-open (PL) Model is easy to pull open and close

	Size in cm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
A	10 x 7	3.94 x 2.76	PL1007
	15 x 7	5.91 x 2.76	PL1507
	20 x 7	7.87 x 2.76	PL2007
	25 x 7	9.84 x 2.76	PL2507

A


Push-to-open (PS) Model is easy to push open and close

	Size in cm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
B	10 x 7	3.94 x 2.76	PS1007
	15 x 7	5.91 x 2.76	PS1507
	20 x 7	7.87 x 2.76	PS2007
	25 x 7	9.84 x 2.76	PS2507

B


U-Shape (US) Model is only with the Side Bars

	Size in cm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
C	10 x 7	3.94 x 2.76	US1007
	15 x 7	5.91 x 2.76	US1507
	20 x 7	7.87 x 2.76	US2007
	25 x 7	9.84 x 2.76	US2507
	25 x 14	9.84 x 5.51	US2514

C



Silicone Mats




**INSTRUMENTS NOT INCLUDED*

The Kent Surgical Silicone Instrument Mat is made from high-grade silicone, offering a durable and reliable solution for instrument protection during the demanding cycles of the Centralized Sterile Supply Department (CSSD). Featuring a unique perforated design, this innovative mat facilitates optimal water and steam circulation, ensuring thorough cleaning and sterilization while minimizing drying time. Additionally, strategically placed raised pins provide universal compatibility for instruments of various dimensions, preventing movement and potential damage during processing.

Size in mm (L x W x H)	Size in Inches (L x W x H)	Catalog Number
520 x 220	20.47 x 8.66	SM5222
460 x 220	18.11 x 8.66	SM4622
380 x 220	14.96 x 8.66	SM3822
260 x 220	10.24 x 8.66	SM2622

SPRI Basket




	Category	Description	Catalog Number
A	Baskets	SPRI Basket Double Frame SS, Stackable and Nestable. L x W x H : 600 x 400 x 100 mm / 23.62 x 15.75 x 3.94 Inch	SP604010
B	Baskets	SPRI Basket Double Frame SS, Stackable and Nestable. L x W x H : 600 x 400 x 150 mm / 23.62 x 15.75 x 5.90 Inch	SP604015
C	Baskets	SPRI Basket Double Frame SS, Stackable and Nestable. L x W x H : 600 x 400 x 200 mm / 23.62 x 15.75 x 7.87 Inch	SP604020

Metal Tray Tags

Kent Surgical offers versatile stainless steel tags that attach securely to the edges of trays or containers. Designed with a flat surface, they accommodate barcodes, labels, or other identification markings. Offered in five sizes, these tags are suitable for steam sterilization and are available either as single units or in packs of 10.




	Size (mm)	Size (inches)	Quantity	Catalog Number
A	35 x 35	1.38 x 1.38	1	MT-3535
			10	MT-3535-10
B	55 x 30	2.17 x 1.18	1	MT-5530
			10	MT-5530-10
C	70 x 30	2.76 x 1.18	1	MT-7030
			10	MT-7030-10
D	90 x 30	3.54 x 1.18	1	MT-9030
			10	MT-9030-10
E	90 x 50	3.54 x 1.97	1	MT-9050
			10	MT-9050-10
F	PERMANENT MARKING ON TAG (PRINT ONLY)			MT-MARK



Peel Pouch Sterilization Rack

These stainless steel racks are engineered to keep peel pouches upright, ensuring optimal contact with the sterilant. They come with slots specifically designed for peel-pouched items, enhancing efficiency in both assembly and drying following steam sterilization



	Category	Description	Catalog Number
A	Peel Pouch Sterilization Rack	Peel Pouch Sterilization Rack 4 Slots L x W : 510 x 150 mm / 20.08 x 5.90 Inch	PSR04
B	Peel Pouch Sterilization Rack	Peel Pouch Sterilization Rack 8 Slots L x W : 230 x 100 mm / 9.06 x 3.94 Inch	PSR05
C	Peel Pouch Sterilization Rack	Peel Pouch Sterilization Rack 5 Slots L x W : 510 x 150 mm / 20.08 x 5.90 Inch	PSR08
D	Peel Pouch Sterilization Rack	Peel Pouch Sterilization Rack 10 Slots L x W : 460 x 100 mm / 18.11 x 3.94 Inch	PSR10
E	Peel Pouch Sterilization Rack	Peel Pouch Sterilization Rack 20 Slots L x W : 410 x 150 mm / 16.14 x 5.90 Inch	PSR21





KENT SURGICAL

Products

- Surgical Instruments
- Surgical Instrument Management Software
- Baskets, Stringers and Accessories

Services

- Repair and Refurbishment of Surgical Instruments



Kent Surgical Industries LLC

Plot No.12, 54th Street, Companies Camp, Al Ain, Abu Dhabi, United Arab Emirates
P: +9713753 2555 | E: office@kentsurgical.com | W: www.kentsurgical.com